


Steps to work with Anea's digital height and weight measuring device

- Placement on the scale
- Carrying out weighing operations
- Height measurement automatically
- Display height, weight and BMI on the screen
- Reporting by scanning the QR code of the display and registering in the Anea Cloud system



Anea Cloud Software



Anea Cloud system

- Unlimited storage of performed analysis
- Show the health status of the person
- Chart of height, weight and BMI changes
- The ability to create a user profile
- Display BMI (body mass index)
- Weight display
- Display BMR (Basal Metabolism)
- Height display
- Body norm charts

Different parts of the software

BMI (body mass index):

- Graph of individual BMI changes in different periods of time
- BMI display and analysis

Display weight:

- Analyzing the individual's weight and displaying the ideal weight
- Chart of weight changes

height display:

- Analysis of a person's height compared to other people of the same gender and age
- Chart of height changes

Show BMR

Growth Chart

